

# BIOVET<sup>®</sup>

# 2023

In partnership with



## Test kits

Detection of allergens, food pathogens, Mycotoxins, GMOs, and other Neogen products



1-888-824-6838

[order@biovet-inc.com](mailto:order@biovet-inc.com)

# Summary

Summary .....	2
Term and conditions .....	3
About Biovet .....	3
Allergens .....	4
Reveal.....	4
Alert.....	5
Veratox.....	6
Mycotoxins.....	7
Reveal.....	7
Veratox.....	8
Others* .....	8
GMOs and ruminants.....	9
Water .....	9
Colitag .....	9
Sanitation .....	10
Accuclean Advanced .....	10
AccuPoint Advanced Next Generation.....	11
Others .....	11
Instruments & material.....	12
Readers .....	12
Pipettes .....	12
Others .....	13
Other products available .....	15
Pathogens (ANSR) .....	15
Pathogens (Reveal) .....	15

MKT-036-EN-V.11 — Revised on 2023-06-14



## Term and conditions

A freight charge of \$25.00 is applicable for orders for equipment under \$500.00.

Orders received before Monday 10 am AM will be delivered in the same week. Orders received past this deadline will be deliverable the following week.

Send your order to: [order@biovet-inc.com](mailto:order@biovet-inc.com)

Effective Date January 4, 2023.

## About Biovet

Biovet is proud to offer you unique, quality and innovative products, in collaboration with its partner Neogen. These are powerful tools allowing you to do your tests in-house, for faster results.

Since 1991, Biovet has been operating state-of-the-art laboratories and carrying out microbiological analyses, detection of pathogens, allergens, mycotoxins and GMOs, as well as in the fields of molecular biology and serology (ELISA).

Our laboratory is accredited by the Standards Council of Canada (SCC) for microbiological testing of foods and edibles (see our accreditation scope on our website).

We are multiplying the actions in order to get closer to you with now 2 DIAGNOSTIC SERVICES LABORATORIES in Quebec: SAINT-HYACINTHE and QUEBEC CITY.

Also, our customers appreciate our HIGH AVAILABILITY, since our laboratories are open 6 days and operate 7 days a week.

In addition, Biovet has an efficient collection service that ensures the transport of samples in Quebec as soon as possible, even in the regions.

For further information, feel free to contact your Sales representative.



## Allergens

Code	Description	Price
TRD-1017	Allergen Extraction Kit	\$67,00
TRD-1010	Allergen Environmental Swabbing Kit (100 swabs)	\$170,00

## Reveal

Neogen's Reveal® and Reveal 3D food allergen detection kits are simple, quick and easy for detecting multiple allergens on work surfaces, ingredients and finished products.



Code	Description	Number of tests	Price	For 4 & +
<b>Reveal 3D</b>				
TRD-1170	Reveal 3D <b>Peanut</b> *	10	\$226,00	\$208,00
TRD-1194	Reveal 3D <b>Almond</b> *	10	\$226,00	\$208,00
TRD-1211	Reveal 3D <b>Crustacea</b>	10	\$226,00	\$208,00
TRD-1186	Reveal 3D <b>Gliadin</b>	20	\$453,00	\$418,00
TRD-1411	Reveal 3D <b>Gluten</b> *†	10	\$262,00	\$239,00
TRD-1221	Reveal 3D <b>Total Milk</b>	10	\$226,00	\$208,00
TRD-1228	Reveal 3D <b>Mustard</b>	10	\$226,00	\$208,00
TRD-1193	Reveal 3D <b>Hazelnut</b> *	10	\$226,00	\$208,00
TRD-1247	Reveal 3D <b>Multi Treenut</b>	10	\$357,00	\$329,00
TRD-1185	Reveal 3D <b>Egg</b> *	10	\$226,00	\$208,00
TRD-1187	Reveal 3D <b>Soy</b> *	10	\$226,00	\$208,00
TRD-1284	Reveal 3D <b>Sesame</b>	10	\$226,00	\$208,00
TRD-1175	Reveal 3D <b>Coconut</b> *	10	\$262,00	\$239,00
Equipment needed but sold separately from the kit				
TRD-1410	* Reveal 3D Food Buffer		\$30,00	
TRD-1409	† Reveal 3D Gluten Food Buffer		\$32,00	
TRD-1408	* Tubes 50 ml graduated polypropylen (25)		\$56,00	

\* Available to analyze on food and ingredients

Code	Description	Number of tests	Price	For 4 & +
	<b>Reveal</b>			
<u>TRD-1058</u>	Reveal <b>Histamine</b>	24	\$310,00	\$285,00

## Alert

Alert® food allergen detection kits are qualitative ELISA-like screening tests that compare up to 5 samples with a control sample.



Code	Description	Price	For 4 & +
TRD-1095	Alert <b>Histamine</b>	\$167,00	\$151,00
TRD-1018	Alert <b>Total Milk</b>	\$418,00	\$375,00
TRD-1013	Alert <b>Egg</b>	\$418,00	\$375,00
TRD-1116	Alert <b>Sulfite</b> (shrimp only)	\$191,00	\$172,00

### Veratox

Veratox® food allergen detection kits are quantitative ELISA-like tests that can simultaneously measure up to 19 samples based on a curve consisting of five controls.



Code	Description	Price	For 4 & +
TRD-1166	Veratox Allergen Starter Kit	UR	
TRD-1016	Veratox VIP for <b>Peanut</b>	\$869,00	\$800,00
TRD-1014	Veratox <b>Peanut</b>	\$869,00	\$800,00
TRD-1097	Veratox <b>Almond</b>	\$869,00	\$800,00
TRD-1215	Veratox <b>Crustacea</b>	\$869,00	\$800,00
TRD-1196	Veratox <b>Gliadin R5</b>	\$869,00	\$800,00
TRD-1197	Gliadin cocktail solution 25ml	\$43,00	
TRD-1205	Gliadin cocktail solution 100ml	\$156,00	
TRD-1198	Gliadin cocktail solution 500ml	\$683,00	
TRD-1027	Veratox <b>Histamine</b>	\$869,00	\$800,00
TRD-1009	Veratox <b>Total Milk</b>	\$869,00	\$800,00
TRD-1252	Veratox <b>Mustard</b>	\$869,00	\$800,00
TRD-1098	Veratox <b>Hazelnut</b>	\$869,00	\$800,00
TRD-1008	Veratox <b>Egg</b>	\$869,00	\$800,00
TRD-1184	Veratox <b>Soy</b>	\$869,00	\$800,00
TRD-1183	Veratox <b>Sesame</b>	\$869,00	\$800,00

# Mycotoxins

## Reveal

The Reveal® and Reveal Q+ test kits are quantitative tests for the analysis of mycotoxins in grains.



Code	Description	Number of tests	Price	For 4 & +
TRD-1300	Reveal Q MAX <b>DON</b>	25	\$418,00	\$379,00
TRD-1281	Reveal Q+ <b>DON</b>	25	\$393,00	\$357,00
TRD-1263	Reveal Q MAX <b>Ochratoxin</b>	25	\$418,00	\$379,00
TRD-1299	Reveal Q+ <b>Ochratoxin</b>	25	\$418,00	\$379,00
TRD-1413	Reveal Q MAX <b>T2/HT-2</b>	25	\$464,00	\$423,00
TRD-1288	Reveal Q+ <b>T2-HT2</b>	25	\$418,00	\$379,00
TRD-1266	Reveal Q MAX <b>Zearalenone</b>	25	\$464,00	\$423,00
TRD-1264	Reveal Q+ <b>Zearalenone</b>	25	\$464,00	\$423,00
TRD-1265	Reveal Q MAX <b>Aflatoxine</b>	25	\$357,00	\$325,00
TRD-1207	<b>Reveal Q+ Aflatoxine</b>	25	\$346,00	\$314,00

## Mycotoxins

### Veratox

Veratox® mycotoxin test kits are quantitative ELISA-like tests that can be assayed simultaneously up to 19 samples based on a curve consisting of five controls.



Code	Description	Price	For 4 & +
TRD-1037	Veratox <b>Ochratoxine</b>	\$464,00	\$423,00
TRD-1260	Veratox <b>Ochratoxine</b> grain	\$464,00	\$423,00
TRD-1029	Veratox <b>T-2/HT2</b>	\$464,00	\$423,00
TRD-1002	Veratox <b>Zearalenone</b>	\$464,00	\$423,00
TRD-1004	Veratox <b>Aflatoxine</b> HS	\$357,00	\$325,00
TRD-1028	Veratox <b>Aflatoxine</b>	\$357,00	\$325,00
TRD-1001	Veratox <b>DON</b> 5/5	\$464,00	\$423,00
TRD-1092	Veratox <b>DON</b> HS	\$464,00	\$423,00
TRD-1109	Veratox <b>DON</b> 2/3	\$464,00	\$423,00
TRD-1038	Veratox <b>Fumonisine</b>	\$464,00	\$423,00
TRD-1128	Veratox <b>Fumonisine</b> 5/10	\$464,00	\$423,00

### Others\*

Code	Description	Price
TRD-1087	Rock-it for mycotoxines	\$2 085,00
TRD-1287	<b>DON</b> Check Sample (MRM), Low, 100g	\$142,00
TRD-1119	<b>DON</b> Check Sample (MRM), Medium, 100 g	\$137,00
TRD-1110	NeoColumn for <b>Aflatoxine</b> (large) - 50 c	\$1 108,00

\* Check Sample are available for other toxins, for more information, contact your Sales representative



## GMOs and ruminants

Code	Description	Number of tests	Price	For 4 & +
TRD-1012	Reveal for CP4 (Roundup Ready)	50	\$273,00	\$257,00
TRD-1034	Reveal for Ruminant in Feed	25	\$391,00	\$359,00
TRD-1195	Reveal for Ruminant in MBM	25	\$391,00	\$359,00

## Water

### Colitag

Code	Description	Price
TRD-1089	Colitag water tes kit P/A – (20)	\$226,00
TRD-1080	Colitag water tes kit P/A – (100)	\$703,00
TRD-1213	Colitag Comparator	\$27,00



## Sanitation

### Accuclean Advanced

AccuClean® Advanced is a simple test to detect protein residues from liquid or solid food products that have been in contact with the surface to be analyzed. This test is completely autonomous, no instruments are required for reading and it requires only minimal training.



Code	Description	Price
TRD-1171	Accuclean Advanced – Detects ATP (100)	\$459.00

## AccuPoint Advanced Next Generation

NEOGEN® New AccuPoint® Advanced Next Generation (NG) is a complete sanitation monitoring system comprising a handheld instrument, innovative samplers, and data manager software, designed to verify your cleaning processes by measuring adenosine triphosphate (ATP). After cleaning your environment, measuring the amount of ATP present on a surface or in liquids allows you to verify your cleaning processes have been effective.



Code	Description	Prix
TRD-1311	Accupoint Advance Next Generation	\$3 618,00
TRD-1242	AccuPoint Advanced Sampler Surface - (100)	\$489,00
TRD-1243	AccuPoint Advanced Sampler Water - (100)	\$489,00
TRD-1244	AccuPoint Advanced Sampler Access - (100)	\$489,00

## Others

To check the availability and price of these products, contact your Sales representative.

Code	Description
9950	AccuPoint Advanced ALP (50 Samplers+ 5 controls)
9971	AccuPoint Advanced ALP (5 controls)

ALP : alkaline phosphatase

## Instruments & material



### Readers

Code	Description	Price
TRD-1042	Neogen Stat-Fax 4700 Microwell Reader	\$3 414,00
TRD-1292	Raptor Reader (Reader + software)	\$4 747,00
TRD-801	Raptor Printer	\$380,00
TRD-803	Raptor printer paper	\$35,00

### Pipettes

Code	Description	Price
TRD-1033	Pipettor 100 $\mu$ L fixed	\$455,00
TRD-1043	Pipettor 50-300 $\mu$ L, 12 Channels	\$1 527,00
TRD-1067	Pipettor 100 $\mu$ L fixed	\$39,00
TRD-1235	Pipettor 100 $\mu$ L fixed for Q+ (no overflow)	\$156,00
TRD-335	Pipette Tips (1000)	\$60,00
TRD-1216	Pipettor 500 $\mu$ L fixed	\$156,00
TRD-1234	Pipettor 1000 $\mu$ L fixed	\$156,00
TRD-1402	Pipettor 400 $\mu$ L fixed	\$208,00
TRD-1088	Pipettor 20-200 $\mu$ L	\$459,00
TRD-1214	Pipette Tips 100 $\mu$ L - 1ml (1000)	\$112,00
TRD-1414	Pipettor 250 $\mu$ L fixed	UR
TRD-1401	Pipettes 400 $\mu$ L exact volume for Raptor (50)	\$19,00
TRD-1236	Sterile Pipettor 10ml (100)	\$63,00

## Others

Code	Description	Price
TRD-1031	Filter Syringes (100)	\$127,00
TRD-1032	Well Holder	\$14,00
TRD-1045	Plastic Wash Bottle, 500 ml (package of 4)	\$37,00
TRD-1047	Incubator with Thermometer	\$738,00
TRD-1049	Microcentrifuge Tube 2 ml (500)	\$105,00
TRD-1056	Whirl-Pak Filter Bags 24oz (250)	\$195,00
TRD-1059	Shaker Water Bath	\$6 582,00
TRD-1066	Agri-Grind	\$60,00
TRD-1068	Whatman Filter # 1 (100)	\$39,00
TRD-1069	Thermometer	\$70,00
TRD-1070	Graduated Plastic Cylinder, 250 ml	\$46,00
TRD-1071	Plastic Wash Bottle, 250 ml	\$10,00
TRD-1072	Nalgene Bottle 1 litre	\$19,00
TRD-1073	Whatman Filter # 4 11cm (100)	\$39,00
TRD-1091	Electronic Digital Scale	\$299,00
TRD-1094	Premoistened Sampling Sponges	\$343,00
TRD-1096	Timer 3-time	\$74,00
TRD-1101	Allergen Wash Buffer Concentrate	\$34,00
TRD-1105	Tubes, 13 x 100 mm (1000)	\$112,00
TRD-1111	Tubes, 12 x 75 mm (1000)	\$72,00
TRD-1135	Graduated Plastic Cylinder, 50 ml	\$15,00
TRD-1136	Membrane Filters, ISO-GRID, 1600-Grid (box of 100)	\$244,00
TRD-1140	ISO-GRID Filtration Unit	\$105,00
TRD-1142	Blender Jar for mixing slides and lid	\$71,00
TRD-1145	Filtration 6-places	\$2 178,00
TRD-1147	Vacuum Tubing Assembly	\$32,00
TRD-1148	ISO-GRID Filtration Unit Clamp	\$105,00
TRD-1152	Blender, Oster	\$117,00
TRD-1154	Reagent Boats (3/pack)	\$45,00
TRD-1164	Tubes w/caps 3.5ml (100)	\$12,00
TRD-1165	Amber glass bottle 500ml	\$16,00
TRD-1190	Filtration ramp	\$1 377,00
TRD-1191	Environmental Sampling Kit	UR
TRD-1202	Pipette Pump 10ml	\$37,00

## Instruments & material / media

---

Code	Description	Price
TRD-1212	Graduated glass Cylinder 50 ml	\$22,00
TRD-1238	Buffer bottle (72/cs)	\$70,00
TRD-1289	Thermal 2 1/4 Paper	\$27,00
TRD-1290	Rack for Reveal Sample Cups	\$32,00
TRD-1400	Raptor Vortex Cartridges (50)	\$46,00
TRD-1025	Fax Printer Paper	\$27,00
TRD-1039	Fax Printer Paper Stat-Fax (10/pk)	\$57,00
TRD-316	Jars and lids for Raptor (50)	UR
TRD-1081	Sample Containers with sodium thiosulfate 100'	\$108,00

## Other products available

To check the availability and price of these products, contact your Sales representative directly.

### Pathogens (ANSR)

Code	Description
9837	ANSR Complete System
9822	ANSR for <b>E. coli</b> - 96 tests
9815	ANSR for <b>E. coli</b> Enrichment Broth, 500g
9815B	ANSR for <b>E. coli</b> Enrichment Broth, 2 kg
9873	ANSR for <b>Listeria</b> Right No- 96 tests
9824	ANSR for <b>Listeria</b> mono - 96 tests
9814	ANSR <b>Listeria</b> Enrichment Broth 1, 500g (environmental sample)
9880	ANSR <b>Listeria</b> LESS PLUS media, 500g (products and environment)
9870	ANSR for <b>Salmonella</b> - 96 tests
	For Enrichment Broth ( <b>Salmonella</b> ), see your Sales representative

### Pathogens (Reveal)

Neogen's Reveal® and Reveal 2.0 test kits contain everything you need to perform the analysis; including the enrichment medium to detect different pathogens.



Code	Description
9750	Reveal <b>E. coli</b> 0157:H7, 8 h, SC - 20 tests
9751	Reveal <b>E. coli</b> 0157:H7, 20 h, SC - 20 test
9714	Reveal 2.0 <b>E. coli</b> 0157:H7 SC – 20 tests (sample of 65g)
9803	Reveal 2.0 for <b>Salmonella</b> , SC avec RV - 20 tests
9805	Reveal 2.2 for <b>Salmonella</b> , SC avec RV et M-broth - 20 tests
	Reveal 2.0, <b>Listeria</b> One step CS for food - 20 tests
	Reveal 2.0, <b>Listeria</b> One step CS for surface - 20 tests
	Reveal 2.0, <b>Listeria</b> Enrichment 48 h Complete System

CS = Complete System

**1-888-824-6838**

option 3

[order@biovet-inc.com](mailto:order@biovet-inc.com)

4375 Beaudry, Saint-Hyacinthe QC J2S 8W2

T. 450 771-7291 • F. 450 771-4158

© Biovet Inc. – May 2023